## Urine CCA (Schisto) ECO Teste - TR.0301C

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- Adriko, M., et al. 2014. Evaluation of circulating cathodic antigen (CCA) urine-cassette assay as a survey tool for Schistosoma mansoni in different transmission settings within Bugiri District, Uganda. *Acta Trop.* (2014).
- Ayele, B., Erko, B., Legesse, M., Hailu, A. and Medhin, G., 2008. Evaluation of circulating cathodic antigen (CCA) strip for diagnosis of urinary schistosomiasis in Hassoba school children, Afar, Ethiopia. *Parasite*. 15: 69-75.
- Colley, D.G., Binder, S., Campbell, C., King, C.H., Tchuem Tchuenté, L.A., N'Goran, E.K., Erko, B., Karanja, D.M., Kabatereine, N.B., van Lieshout, L., and Rathbun, S. 2013. A Five-Country Evaluation of a Point-of-Care Circulating Cathodic Antigen Urine Assay for the Prevalence of Schistosoma mansoni. *Am. J. Trop. Med. Hyg.*, 88(3), 2013, pp. 426-432.
- Coulibaly, J.T., N'Gesso, Y.K., Knopp, S., N'Guessan, N.A., Silué, K.D., van Dam G.J., N'Goran, E.K. and Utzinger, J. 2013. Accuracy of Urine Circulating Cathodic Antigen Test for the Diagnosis of Schistosoma mansoni in Preschool-Aged Children before and after Treatment. *PLoS Negl Trop Dis* 7(3): e2109. doi:10.1371/journal.pntd.0002109.
- Legesse, M. and Erko, B., 2008. Field-based evaluation of a reagent strip test for diagnosis of Schistosomiasis Mansonii by detecting circulating cathodic antigen (CCA) in urine in low endemic area in Ethiopia. *Parasite*, 15: 151-155.
- Legesse, M. and Erko, B., 2007. Field-based evaluation of a reagent strip test for diagnosis of Schistosoma mansoni by detecting circulating cathodic antigen in urine before and after chemotherapy. *Trans. R. Soc. Trop. Med. Hyg.* 101:668-673.
- Midzi, N., Butterworth, A.E., Muluza, T., Munyati, S., Deelder, A.M., Van Dam, G.J., 2008. The use of circulating cathodic antigen (CCA) strips for the diagnosis of urinary schistosomiasis. *Trans R Soc Trop Med Hyg* (2008), doi:10.1016/j.trstmh.2008.08.018.
- Odogwu S.E., Ramamurthy N.K., Kabatereine N.B., Kazibwe F., Tukaheba E., Webster J.P., Fenwick A., Stothard J.R. 2006. Schistosoma mansoni in infants (aged < 3 years) along the Ugandan shoreline of Lake Victoria. *Ann Trop Med Parasitol*, 100:315-326.
- Pauline N. M. Mwinzi, Nupur Kittur, Elizabeth Ochola, Philip J. Cooper, Carl H. Campbell Jr., Charles H. King and Daniel G. Colley. 2015. Additional evaluation of the point-of-contact circulating cathodic antigen assay for Schistosoma mansoni infection. *Frontiers in Public Health* 2015.
- Shane H.L., Verani J.R., Abudho B., Montgomery S.P., Blackstock A.J., et al., 2011. Evaluation of Urine CCA Assays for Detection of Schistosoma mansoni Infection in Western Kenya. *PLoS Negl Trop Dis* 5(1): e951. doi:10.1371/journal.pntd.0000951.
- Sousa-Figueiredo J.C., Pleasant J., Day M., Betson M., Rollinson D., Montresor A., Kazibwe F., Kabatereine N.B., Stothard J.R. 2010. Treatment of intestinal schistosomiasis in Ugandan preschool children: best diagnosis, treatment efficacy and side-effects, and an extended praziquantel dosing pole. *International Health*, 2:103 - 113.
- Sousa-Figueiredo, J.C., Betson, M., Kabatereine, N.B., Stothard, J.R. 2013. The Urine Circulating Cathodic Antigen (CCA) Dipstick: A Valid Substitute for Microscopy for Mapping and Point-Of-Care Diagnosis of Intestinal Schistosomiasis. *PLoS Negl Trop Dis* 7(1): e2008. doi:10.1371/journal.pntd.0002008.
- Standley C.J., Lwambo N.J.S., Lange C.N., Kariuki H.C., Adriko M., Stothard J.R. 2010. Performance of circulating cathodic antigen (CCA) urine-dipsticks for rapid detection of intestinal schistosomiasis in schoolchildren from shoreline communities of Lake Victoria. *Parasites & Vectors*, 3:7.
- Stothard, J.R., 2009. Improving control of African schistosomiasis: towards effective use of rapid diagnostic tests within an appropriate disease surveillance model. *Trans R Soc Trop Med Hyg*, 103:325 - 332.

- Stothard, J.R., Kabatereine, N.B., Tukaheba, E.M., Kazibwe, F., Rollinson, D., Mathieson, W., Webster, J.P. and Fenwick, A., 2006. Use of circulating cathodic antigen (CCA) dipsticks for detection of intestinal and urinary schistosomiasis. *Acta Trop.* 97:219-228.
- Stothard J.R., Sousa-Figueiredo J.C., Betson M., Adriko M., Arinaitwe M., et al., 2011. Schistosoma mansoni Infections in Young Children: When Are Schistosome Antigens in Urine, Eggs in Stool and Antibodies to Eggs First Detectable? *PLoS Negl Trop Dis* 5(1): e9388. doi:10.1371/journal.pntd.0000938.
- Tchuem Tchuenté, L.-A., Kueté Fouodo, C.J., Kamwa Ngassam, R.I., Sumo, L., Dongmo Noumedem, C., Mérémi Kenfack C., Feussom Gipwe, N., Dankoni Nana, E., Stothard, J.R. and Rollinson, D. 2012. Evaluation of Circulating Cathodic Antigen (CCA) Urine-Tests for Diagnosis of Schistosoma mansoni Infection in Cameroon. *PLoS Negl Trop Dis* 6(7): e1758. doi:10.1371/journal.pntd.0001758.
- Van Dam, G.J., Wichers, J. H., Falcao Ferreira, T. M., Ghati, D., van Amerongen, A. and Deelder, A. M., 2004. Diagnosis of Schistosomiasis by Reagent Strip Test for Detection of Circulating Cathodic Antigen. *J. Clin. Microbiol.* 42, No. 12, Dec. 2004, p. 5458-5461.
- Van Lieshout, L., De Jonge, N., Mansour, M.M., Bassily, S., Kriger, F.W., Deelder, A.M., 1993. Circulating cathodic antigen levels in serum and urine of schistosomiasis patients before and after chemotherapy with praziquantel. *Trans R Soc Trop Med Hyg.* 87:311-2.

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- Coulibaly, J.T., N'Gesso, Y.K., Knopp, S., N'Guessan, N.A., Silué, K.D., et al. 2013. Accuracy of Urine Circulating Cathodic Antigen Test for the Diagnosis of Schistosoma mansoni in Preschool-Aged Children before and after Treatment. *PLoS Negl Trop Dis* 7(3): e2109. doi:10.1371/journal.pntd.0002109.
- Koukounari, A., Donnelly, C.A., Moustaki, I., Tukaheba, E.M., Kabatereine, N.B., et al. 2013. A Latent Markov Modelling Approach to the Evaluation of Circulating Cathodic Antigen Strips for Schistosomiasis Diagnosis Pre- and Post-Praziquantel Treatment in Uganda. *PLoS Comput Biol* 9(12): e1003402. doi:10.1371/journal.pcbi.1003402.
- Lamberton, P.H.L., Kabatereine, N.B., Oguttu, D.W., Fenwick, A., Webster, J.P. 2014. Sensitivity and Specificity of Multiple Kato-Katz Thick Smears and a Circulating Cathodic Antigen Test for Schistosoma mansoni Diagnosis Pre- and Post-repeated-Praziquantel Treatment. *PLoS Negl Trop Dis* 8(9): e3139. doi:10.1371/journal.pntd.0003139.
- Coulibaly, J.T., N'Gesso, Y.K., Knopp, S., Keiser, J., N'Goran, E.K., Utzinger, J. 2012. Efficacy and Safety of Praziquantel in Preschool-Aged Children in an Area Co-Endemic for Schistosoma mansoni and S. haematobium. *PLoS Negl Trop Dis* 6(12): e1917. doi:10.1371/journal.pntd.0001917.
- De Clercq, D., Sacko, M., Vercurysse, J., Vanden Bussche, V., Landouré, A., Diarra, A., Grysels, B., Deelder, A. 1997. Assessment of cure by detection of circulating antigens in serum and urine, following schistosomiasis mass treatment in two villages of the Office du Niger, Mali. *Acta Trop.* 68:339-46.
- Kremsner, P.G., Enyong, P., Kriger, F.W., De Jonge, N., Zotter, G.M., Thalhammer, F., Mühlischlegel, F., Bienzie, U., Feldmeier, H., Deelder, A.M. 1994. Circulating anodic and cathodic antigen in serum and urine from Schistosoma haematobium-infected Cameroonian children receiving praziquantel: a longitudinal study. *Clin Infect Dis*, 18:408-13.
- Van Lieshout, L., De Jonge, N., Mansour, M.M., Bassily, S., Kriger, F.W., Deelder, A.M. 1993. Circulating cathodic antigen levels in serum and urine of schistosomiasis patients before and after chemotherapy with praziquantel. *Trans R Soc Trop Med Hyg.* 87:311-2.

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